


AddOn[®] hydrophilic · toric-trifocal



Technical specifications

Refractive Index 1.46

Abbe No. 58

Lens type One piece sulcus lens

Material Biocompatible hydrophilic acrylic (25%)

Haptic 4 Flex-Haptic

Optic spheric

Optic form convex-concave


Filter UV filter

Sterilisation Steam

Injection system 1stInjct 2.0HB, 1stInjct 2.4HB

■ Basis function

Optical Power (spherical equivalent)




■ Cylinder function

Cylindrical Power



■ Trifocal function

Addition:



AddOn[®] hydrophilic · toric-trifocal

Article No.	SEQ	Steps	Cylinder	Steps	Addition
A4HW0M	-3.0 D to +3.0 D	0.5 D	1.00 D to 4.50 D	0.50 D	+3.00 D



For the latest version,
scan the QR code.
1stq.de/addon1053

**Do you have any
questions?**
info@1stq.de



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en