


AddOn[®] SML hydrophilic



Technical specifications

Refractive Index	1.46
------------------	------

Abbe No.	58
----------	----

Lens type	One piece sulcus lens
-----------	-----------------------

Material	Biocompatible hydrophilic acrylic (25%)
----------	---

Haptic	4 Flex-Haptic
--------	---------------

Optic	spheric
-------	---------

Optic form	convex-concave
------------	----------------

Filter	UV filter
--------	-----------

Sterilisation	Steam
---------------	-------

Injection system	1stInject 2.0HB, 1stInject 2.4HB
------------------	----------------------------------

Magnifier function

Central addition:

+10.0 D

AddOn[®] SML hydrophilic

Article No.	Powers	Addition
A4MW00	On request	+10.00 D

AddOn[®] SML (Scharioth Macula lens)

AddOn[®] SML is an effective treatment option for patients with advanced maculopathy. The central optic of +10 D addition at 1.5 mm diameter acts as a miniaturised magnifying glass.

Why AddOn[®] for AMD?

The SML uses the special advantages of the AddOn[®]:

- Standardised implantation
- Independence from the lens status of the pseudophakic eye
- Distance visual acuity and visual field remain unchanged
- Reversibility in case of intolerance
- Significantly reduced stigmatisation of the AMD patient

Indication

- Dry AMD from a visual acuity of 0.1 (monocular implantation in the dominant eye)
- Myopic macular degeneration
- Diabetic macular degeneration and other maculopathies

Contraindication

- Wet AMD
- Severe zonulopathy
- Subluxation
- Aphakia
- Progressive glaucoma
- Iris neovascularisation
- Poor retinal function



For the latest version, scan the QR code.

[1stq.de/addon1050](https://www.1stq.de/addon1050)

Do you have any questions?

info@1stq.de



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany

www.1stq.de/en