

Technical specifications

Refractive Index	1.46
------------------	------

Abbe No.	58
----------	----

Lens type	One piece sulcus lens
-----------	-----------------------

Material	Biocompatible hydrophilic acrylic (25%)
----------	---

Haptic	4 Flex-Haptic
--------	---------------

Optic	spheric
-------	---------

Optic form	convex-concave
------------	----------------


Filter	UV filter
--------	-----------

Sterilisation	Steam
---------------	-------

Injection system	1stInject 2.0HB, 1stInject 2.4HB
------------------	----------------------------------

■ Basis function

Optical Power



■ EDOF function

low addition

AddOn[®] hydrophilic · EDOF

Article No.	Powers	Steps	Addition
A4DW0H	0.0 D		low addition
A4EW0H	-5.0 D to -0.5 D	0.25 D	low addition
A4EW0H	+0.5 D to +5.0 D	0.25 D	low addition



For the latest version,
scan the QR code.
1stq.de/addon1069

**Do you have any
questions?**
info@1stq.de



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en