


AddOn[®] hydrophilic · toric



Technical specifications

Refractive Index	1.46
------------------	------

Abbe No.	58
----------	----

Lens type	One piece sulcus lens
-----------	-----------------------

Material	Biocompatible hydrophilic acrylic (25%)
----------	---

Haptic	4 Flex-Haptic
--------	---------------

Optic	spheric
-------	---------

Optic form	convex-concave
------------	----------------

Filter	UV filter
--------	-----------

Sterilisation	Steam
---------------	-------

Injection system	1stInjct 2.0HB to Cyl. 4.50 D 1stInjct 2.4HB
------------------	---


■ Basis function

Optical Power (spherical equivalent)



■ Cylinder function

Cylindrical Power



AddOn[®] hydrophilic · toric

Article No.	SEQ	Steps	Cylinder	Steps
A4TW0T	0.0 D		1.00 D	0.50 D
A4TW0T	0.0 D		1.50 D to 4.50 D	0.75 D
A4TW00	0.0 D		5.25 D to 8.25 D	0.75 D
A4TW00	0.0 D		9.00 D to 11.00 D	1.00 D
A4FW0T	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	1.00 D	0.50 D
A4FW0T	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	1.50 D to 4.50 D	0.75 D
A4FW00	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	5.25 D to 8.25 D	0.75 D
A4FW00	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	9.00 D to 10.00 D	1.00 D
A4FW00	-3.0 D to +8.0 D	0.25 D	11.00 D	



For the latest version,
scan the QR code.
[1stq.de/addon1052](https://www.1stq.de/addon1052)

**Do you have any
questions?**
info@1stq.de



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en